

# INTRAVENOUS (IV) REMINDER



- Advanced beeper warns when intravenous (IV) bag is empty or jammed to prevent backflow of blood and clotting.
- I.R sensor design, reliable way to monitor IV therapy.
- Fully compatible IV drip chamber sizes, easy to install and remove.
- Cable-less/compact design easy to carry.



- Battery life- over 500 hours.
- L.E.D status indicator for low battery.
- Beeper warns when IV drip is jammed or empty.
- Wireless or cable monitor options available.

Clip concept design, easy to install.



If tubing is jammed or IV bag is empty, alarm beeper will warn and show RED coloured light.



When IV bag is empty.



When tubing is jammed.