

LIGHT WHERE YOU NEED IT



**14K**  
LUX AT 4"



**3 HR.**  
RUN TIME



**IPX5**  
RESISTANT



**MADE IN  
THE USA**

**STERILE EO**  
STERILIZED USING  
ETHYLENE OXIDE



**NOT MADE  
WITH  
NATURAL  
RUBBER  
LATEX**



**INDIVIDUALLY  
PACKAGED**



**SINGLE USE  
DISPOSABLE**

**Scintillant® Suction**



**Part No. 2658-01-0002**

The Scintillant® Suction is a handheld surgical light assembled to a Yankauer style suction tip designed to combine fluid aspiration with bright LED illumination. The Suction is ready for use directly from the package, enabling users to immediately navigate obstructions safely, clearly visualize anatomy, and aspirate fluids easily.

CE 0344 3-Year Shelf Life

**Scintillant® Retractor**



**Part No. 2658-01-0004**

The Scintillant® Retractor is a handheld surgical light affixed to a low profile flat retractor designed to combine moderate blunt retraction of tissue with bright LED illumination. The ergonomic design enables users to safely navigate critical obstructions, clearly visualize anatomy and efficiently retract tissue. The Retractor is ready for use immediately upon removal from its protective packaging.

CE 0344 3-Year Shelf Life

**Scintillant® Zip-Strap**



**Part No. 2658-01-0001**

The Scintillant® Zip-Strap is a handheld surgical light featuring rubber straps for simple attachment to surgical instruments. The non-tethered construction enables users to enhance operating room instruments, clearly illuminate the surgical field and safely navigate critical obstructions.

CE 0344 3-Year Shelf Life

# Scintillant® SURGICAL LIGHT

<b>SURGICAL PROCEDURE</b>	<b>Suction</b>	<b>Retractor</b>	<b>Zip-Strap</b>
<b>GYNECOLOGICAL</b>	<b>2658-01-0002</b>	<b>2658-01-0004</b>	<b>2658-01-0001</b>
Abdominal Hysterectomy	●		●
Abdominal Hysterectomy-Obese Pt	●	●	●
Pelvic Floor Reconstruction	●	●	●
Prolapse Procedures	●	●	●
Vaginal Hysterectomy	●	●	●
Vaginal Vault Suspension	●	●	●
<b>UROLOGIC</b>			
Bladder Suspension	●	●	
Cystectomy	●		●
Nephrectomy	●		●
Prostatectomy	●	●	●
Renal Decapsulation			●
<b>COLORECTAL</b>			
Abdominal Perineal Resection	●		●
Lower Anterior Resection	●	●	●
Sigmoid Colon Resections	●	●	●
<b>GENERAL</b>			
Hepatic Procedures	●	●	●
Lymph Node Dissection	●	●	●
Oncological Procedures	●	●	●
<b>ORTHOPEDIC</b>			
Hip Revision	●		●
Total Hip Replacement	●		●
<b>PLASTIC/RECONSTRUCTIVE</b>			
Abdominoplasty	●	●	●
Breast Augmentation	●	●	●
Craniofacial Reconstruction	●	●	●
Reconstruction Mammoplasty	●	●	●
Thoracic	●	●	●
Esophagectomy	●	●	●
Lobectomy	●	●	●
Nipple Sparing Mastectomy	●	●	●
Thoracotomy	●	●	●
<b>SPINE</b>			
Deep Cavity Illumination			●

SEE IFU FOR INTENDED USE AND CONTRAINDICATIONS

2658-80-0044 Rev A